

מחנה
העבר



Heat Exchangers in Pharmaceutical Water Systems

- Guidance, Challenges & Solutions
8th November 2021

Per-Åke Olsson

Global Industry Manager
Alfa Laval

Amir Jahangiri

Global Product Manager
Alfa Laval

Per-Åke Olsson

- Introduction of speaker



Per-Åke Olsson
Global Industry Manager, Pharma & Biotech, Alfa Laval

Alfa Laval AB, 2002- (Lund, Shanghai, Lund)
Astra Zeneca R&D Lund, Sweden, 1997-2002
NiMe Hydrid AB, Mönsterås, Sweden, 1993-1997

Member of the Pharmaceutical Technology Europe editorial advisory board between 2006 and 2008, and has been a speaker and chairman at several BioPharma conferences and seminars in Asia, America and Europe

MSc, Mechanical Engineering, University of Lund, Sweden 1987-1993
eMBA, University of Warwick, UK 2000-2004

telephone: +46 722 171 226

email: perake.olsson@alfalaval.com.

Amir Jahangiri

- Introduction of speaker



Amir Jahangiri

Global Product Mgt, BU FHT

Pharmaceutical HEX; Pharma-line and Pharma-X , Pharma-line Point of Use , Alfa Laval

Alfa Laval AB, 2018- (Lund)

PolyPeptide Group , Sweden, 2012-2017

Novozymes , Sweden, 2011-2012

Experienced in Lean production and six sigma in pharma industry. Responsible for development of heat exchangers for WFI /PW application.

BSc, Chemical Engineering, Petroleum University , Iran 1991-1996

MBA, Industrial Management Institute , Iran 2004-2006

Food and Drug program , Folk University , Sweden 2010-2011

telephone: +46 722 231 087

email: amir.jahangiri@alfalaval.com.

Patrik Arvidsson

- Introduction of speaker



Patrik Arvidsson
Senior Specialist Pharma
Industry responsible heat exchangers, Biotech & Pharmaceutical, Alfa Laval

Alfa Laval, 1998-

Responsible for development of heat exchangers to Pharma. Expert in heat transfer and pharma demands, and has been a speaker at several BioPharm conferences and seminars in Asia, America and Europe.

MSc, Chemical Engineering, Food and pharma engineering, University of Lund, Sweden 1990-1995
Tech. Lic. University of Lund, Sweden 1995-1998.

patrik.arvidsson@alfalaval.com

Agenda

– What will be covered



Part 1:

- Introduction to water storage & distribution loops
- What are the challenges and what does the guidance say about heat exchanger for water storage & distribution loop?
 - FDA guide to inspection
 - ISPE
 - ASME BPE

Part 2:

- How Pharma-line S&T and Pharma-line Point of Use fulfil the guidance and solve the challenges

USP water specification



	WFI	PW
Conductivity $\mu\text{S}/\text{cm}$ @ 25 °C*	1,3	1,3
TOC, ppb	< 500	< 500
Microbial, cfu/100ml	10	10 000**
Endotoxin, max	0,25 EU/ml	Not spec.

* Conductivity levels depending on measured temperature

** Purified water specified as 100 cfu/ml

Water For Injection Monographs					
Attribute	USP	EP	JP	ChP	IP
Source water	US, EU, Japan, WHO drinking water	Human consumption	JP water specification	Potable water or Purified water	Potable water or Purified water
Production Method	Distillation or suitable processes	Distillation	Distillation or RO with UF from Purified Water	Distillation	Distillation
Microbial (cfu/100 ml)	10	10	10	10	10
Conductivity ($\mu\text{S}/\text{cm}$ at 25°C)	1,3 (3 stage)	1,3 (3 stage)	1,3 on-line or 2,1 off-line	1,3 (3 stage)	1,3 (3 stage)
TOC (mg/l)	0,5	0,5	0,5 (0,3 for control)	0,5	0,5
Endotoxins (EU/ml)	0,25	0,25	0,25	0,25	0,25
Nitrates (ppm)	-	0,2	-	0,2	Required
Acidity/Alkalinity	-	-	-	Required	-
Ammonium (ppm)	-	-	-	Required	-
Oxidizable substances	-	Required	-	0,2	-


Purified Water Monographs					
Attribute	USP	EP	JP	ChP	IP
Source water	US, EU, Japan, WHO drinking water	Human consumption	JP water specification	Potable water or Purified water	Potable water or Purified water
Production Method	Suitable processes	Suitable processes	Distillation, ion-exchange, UF, or combination	Distillation, ion-exchange, or suitable processes	Distillation, ion-exchange, or suitable processes
Microbial (cfu/ml)	100	100	100	100	100
Conductivity ($\mu\text{S}/\text{cm}$ at 25°C)	1,3 (3 stage)	5,1 (1 stage)	1,3 on-line or 2,1 off-line	5,1 (1 stage)	1,3 (3 stage)
TOC (mg/l)	0,5	0,5 (optional)	0,5 (0,3 for control)	0,5	0,5
Nitrates (ppm)	-	0,2	-	0,2	Required
Acidity/Alkalinity	-	-	-	-	Required
Ammonium (ppm)	-	-	-	0,2	Required
Oxidizable substances	-	Required	-	-	-

Storage & distribution loop

Water Generation




Storage & Distribution



Ambient or cold water systems

Hot water systems



Challenges & Guidance's for Heat exchangers



- Pharmaceutical water systems

Challenges:

1. Cross contamination
2. Sanitization
3. Velocity
4. Surface finish
5. Welding
6. Passivation
7. Rouging
8. Installation & operation
9. Service & maintenance

Guidance's:

- US FDA, Inspection guides
<https://www.fda.gov/inspections-compliance-enforcement-and-criminal-investigations/inspection-references/inspection-guides>
- ASME BPE – 2019
- ISPE Baseline Guide Volume 4,
Water and Steam Systems, Third edition, 2019



1. Cross contamination

- ISPE & FDA

Two ways to control:

- Over pressure (Plate or single tube sheet heat exchangers)
- Double tube sheet S&T or Tube-in-tube heat exchangers

Recommendations:

- Hygienic applications:
 - Double tube sheet S&T or Tube-in-tube
- Pretreatment or non-hygienic:
 - Plate or single tube sheet S&T

ISPE Vol 4

- 8.3.4.2 Heat Exchangers – Design Considerations

US FDA Inspection Guides

- GUIDE TO INSPECTIONS OF HIGH PURITY WATER SYSTEMS
- HEAT EXCHANGERS TO AVOID CONTAMINATION

2. Sanitization & Bacteria control

- ISPE & FDA



It is recognized that hot (65 - 80°C) systems are self sanitizing

Sanitization:

- S&T, no need for special arrangement
- Gasketed plate, heat sanitization ok, chemical sanitization requires additional rinse

Drainability:


- S&T heat exchangers, Drainability with slope or for U-tubes with weep holes at the low point
- Plate heat exchangers, not drainable (Gasketed can be disassembled for drainability)



non-recirculating water systems be drained daily and water not be allowed to sit in the system.

2. Sanitization & Bacteria control

- Bacteria found in pharmaceutical water systems



Bacterium	Optimum (°C) /Upper (°C)
Escherichia coli	37/45
Staphylococcus aureus	30-37/45
Pseudomonas maltophilia	35/41
Pseudomonas aeruginosa	37/42
Pseudomonas fluorescens	25-30/39
Listeria monocytogenes	30-37/45
Campylobacter jejuni	37-42/45
Clostridium perfringens	37/50
Shigella spp.	35-37/47

Belong to Mesophile class

Growth rate vs temperature for five environmental classes of prokaryotes


1 & 2. ASME BPE - SD-3.6 Heat Exchanger Equipment

- Design for Cross contamination and Sanitization

- Easy to clean
- Inspectable
- Bending process
- Sample
- Free of liquid penetrant indications
- Double tube sheet (Figure SD- 3.6.1-1)
- Orientation/Drainability
- Thermal and Mechanical calculations (Operation & SIP)
- Design pressure: Process side \geq Utility side
- Cleaning & Steaming
- Gaskets and Seals

(f) Minimum recommended bend radii for heat exchangers should be as follows:

Nominal Tube O.D.		Minimum Bend Radius	
in.	mm	in.	mm
0.375	9.5	0.625	15.2
0.500	12.7	0.750	19.1
0.625	15.8	0.938	23.8
0.750	19.1	1.125	28.6
1.000	25.4	1.500	38.1




3. Velocity


- ISPE



- Turbulent flow, Reynolds > 4000 (Biofilm)



Re $< 4\ 000$ – Laminar flow
Low shear forces



↳ Boundary layer

Re $> 4\ 000$ – Turbulent flow
Higher shear forces

Velocity	Boundary layer
0,9 m/s (3 ft/s)	120 μm
2,4 m/s (8 ft/s)	50 μm
Very high velocity	10-20 μm

4. Surface finish

- ISPE & ASME BPE

- Ra 0,3-0,8 μm (10-30 μinch)
- Mechanically or Electro polished
- Polishing rests
- Passivation

Table SF-2.4.1-1 R_a Readings for Metallic Process Contact Surfaces

Surface Designation	Mechanically Polished [Note (1)]	
	R_a Max.	
	$\mu\text{in.}$	μm
SF0	No finish requirement	No finish requirement
SF1	20	0.51
SF2	25	0.64
SF3	30	0.76
Electropolished		
R_a Max.		
	$\mu\text{in.}$	μm
SF4	15	0.38
SF5	20	0.51
SF6	25	0.64

GENERAL NOTES:

- All R_a readings are to be in accordance with ASME B46.1.
- All R_a readings are taken across the lay, wherever possible.
- No single R_a reading shall exceed the R_a max. value in this table.
- Other R_a readings are available if agreed on between the owner/user and supplier, not to exceed values in this table.

NOTE:

- Or any other finishing method that meets the R_a max.

Table SF-2.2-1 Acceptance Criteria for Metallic Process Contact Surface Finishes

Anomaly or Indication	Acceptance Criteria
Pits/porosity	If diameter <0.020 in. (0.51 mm) and bottom is shiny [Notes (1) and (2)]. Pits <0.003 in. (0.08 mm) diameter are irrelevant and acceptable.
Cluster of pits/porosity	No more than 4 pits per 0.5 in. (13 mm) \times 0.5 in. (13 mm) inspection window. The cumulative total diameter of all relevant pits shall not exceed 0.040 in. (1.02 mm).
Dents	None accepted [Note (3)]
Finishing marks	If R_a max. is met
Welds	Welds used in the as-welded condition shall meet the requirements of MJ-8. Welds finished after welding shall be flush with the base metal, and concavity and convexity shall meet the requirements of MJ-8. Such finishing shall meet the R_a requirements of Table SF-2.4.1-1.
Nicks	None accepted
Scratches	For tubing, if cumulative length is <12.0 in. (305 mm) per 20 ft (6.1 m) tube length or prorated and if depth is <0.003 in. (0.08 mm) For fittings, valves, and other process components, if cumulative length is <0.25 in. (6.4 mm), depth <0.003 in. (0.08 mm), and R_a max. is met For vessels, if length <0.50 in. (13 mm) at 0.003 in. (0.08 mm) depth and if <3 per inspection window [Note (4)]
Surface cracks	None accepted
Surface inclusions	If R_a max. is met
Surface residuals	None accepted, visual inspection
Surface roughness (R_a)	See Table SF-2.4.1-1
Weld slag	For tubing, up to 3 per 20 ft (6.1 m) length or prorated, if <75% of the width of the weld bead For fittings, valves, vessels, and other process components, none accepted (as welded shall meet the requirements of MJ-8 and Table MJ-8.4-1)
Blistering	None accepted

GENERAL NOTE: This table covers surface finishes that are mechanically polished or any other finishing method that meets the R_a max.

NOTES:

- Black bottom pit of any depth is not acceptable.
- Pits in superaustenitic and nickel alloys may exceed this value. Acceptance criteria for pit size shall be established by agreement between owner/user and supplier. All other pit criteria remain the same.
- For vessels, dents in the area covered by and resulting from welding dimple heat transfer jackets are acceptable.
- An inspection window is defined as an area 4 in. \times 4 in. (100 mm \times 100 mm).

Table SF-2.2-2 Additional Acceptance Criteria for Electropolished Metallic Process Contact Surface Finishes

Anomaly or Indication	Acceptance Criteria
Cloudiness	Acceptable if R_a max. is met
End grain effect	Acceptable if R_a max. is met
Fixture marks	Acceptable if electropolished
Haze	Acceptable if R_a max. is met
Interrupted electropolish	Acceptable if R_a max. is met
Orange peel	Acceptable if R_a max. is met
Stringer indication	Acceptable if R_a max. is met
Weld whitening	Acceptable if R_a max. is met
Variance in luster	Acceptable if R_a max. is met

5. Welding

- ASME BPE, Chapter 5, Part MJ Materials joining

MJ-1 Purpose & Scope

MJ-2 Materials

MJ-2 Joint design and preparation

MJ-4 Procedure qualifications

MJ-6 Performance qualifications

MJ-7 Examination, inspection, and testing

MJ-8 Acceptance criteria

MJ-9 Joining of polymeric materials

MJ-10 Documentation requirements

MJ-11 Passivation




Figure 6. Flow chart for orbital welding/inspection/documentation of stainless steel welds.

5. Welding

– ASME BPE MJ-8 Acceptance criteria

- Cracks
- Lack of fusion
- Incomplete penetration
- Porosity
- Inclusions
- Undercut
- Groove weld concavity
- Fillet weld concavity
- Discoloration (HAZ)
- Discoloration (weld bead)

- Oxide island
- Reinforcement
- Tack welds
- Arc strikes
- Overlap
- Weld bead width
- Minimum fillet weld size
- Misalignment

Table MJ-8.2-1 Visual Examination Acceptance Criteria for Welds on Metallic Pressure Vessels and Tanks

Discontinuities	Welds on Process Contact Surfaces			Welds on Non-Process Contact Surfaces	
	Welds Left in the As-Welded Condition	Prior to Postweld Finishing	After Postweld Finishing	Welds Left in the As-Welded Condition	After Postweld Finishing
Cracks	None	None	None	None	None
Lack of fusion	None	None	None	None	None
Incomplete penetration	None on process contact side; otherwise, see Note (1)	None on process contact side; otherwise, see Note (1)	None on process contact side; otherwise, see Note (1)	See Notes (1) and (2)	See Notes (1) and (2)
Porosity	None open to the surface; otherwise, see Note (1)	See Note (1)	See Table SF-2.2-1 for acceptance criteria for pits/porosity	None open to the surface; otherwise, see Note (1)	None open to the surface; otherwise, see Note (1)
Inclusions [metallic (e.g., tungsten) or nonmetallic]	None open to the surface; otherwise, see Note (1)	See Note (1)	None open to the surface; otherwise, see Note (1)	None open to the surface; otherwise, see Note (1)	None open to the surface; otherwise, see Note (1)
Undercut	None	See Note (1)	None	See Note (1)	See Note (1)
Groove weld concavity	See Note (1)	See Note (1)	Maximum of 10% of the nominal wall thickness of thinner member	See Note (1)	See Note (1)
Fillet weld convexity	1/16 in. (1.5 mm) max.	Per applicable design and fabrication code	1/32 in. (0.8 mm) max.	See Note (1)	See Note (1)

6. Passivation

- ISPE & ASME BPE

- Re-passivation (ISPE)
- Passivated before being placed in service (ASME BPE)
- Passivation of electropolished surfaces is not required unless the surface has been altered (ASME BPE)
- ASME BPE Nonmandatory Appendix E, Passivation procedure qualification

NONMANDATORY APPENDIX E PASSIVATION PROCEDURE QUALIFICATION		Table E-2.2-1 Minimum Surface Treatments for Process Qualification Schemas		Table E-2.2-2 Passivation Process		Table E-2.2-3 Passivation Process (E-2)		Table E-2.2-4 Test Matrix for Evaluation of Cleaned and/or Passivated Surfaces		Table E-2.2-5 Test Matrix for Evaluation of Cleaned and/or Passivated Surfaces (E-2)		Table E-2.2-6 Test Matrix for Evaluation of Cleaned and/or Passivated Surfaces (E-2)		Table E-2.2-7 Test Matrix for Evaluation of Cleaned and/or Passivated Surfaces (E-2)							
Process Type	Process Description	Component Size (L)	Component Size (D)	Component Size (R)	Component Size (T)	Process Type	Process Description	Type of Test	Test Description	Type of Test	Test Description	Type of Test	Test Description	Type of Test	Test Description						
E-1 GENERAL	This Appendix provides basic information and offers guidelines for surface passivation, equipment qualification, and process control. It is intended to be used in conjunction with the requirements of 3.1.1.1 and 3.1.1.2. The minimum surface treatment and passivation procedures are intended to be used in conjunction with the requirements of 3.1.1.1 and 3.1.1.2. The minimum surface treatment and passivation procedures are intended to be used in conjunction with the requirements of 3.1.1.1 and 3.1.1.2.	<p>E-1.1 Passivation Procedure (E-1)</p> <p>The passivation procedure shall utilize wettable and non-wettable surfaces. The passivation procedure shall be performed in accordance with the requirements of 3.1.1.1 and 3.1.1.2. The passivation procedure shall be performed in accordance with the requirements of 3.1.1.1 and 3.1.1.2.</p>				<p>E-1.2 Certificate of Passivation Compliance</p> <p>The passivation procedure shall comply with the requirements of 3.1.1.1 and 3.1.1.2. The passivation procedure shall be performed in accordance with the requirements of 3.1.1.1 and 3.1.1.2. The passivation procedure shall be performed in accordance with the requirements of 3.1.1.1 and 3.1.1.2.</p>				<p>E-1.3 Acceptance Criteria for Cleaned and/or Passivated Surfaces (E-2)</p> <p>The passivation procedure shall be performed in accordance with the requirements of 3.1.1.1 and 3.1.1.2. The passivation procedure shall be performed in accordance with the requirements of 3.1.1.1 and 3.1.1.2. The passivation procedure shall be performed in accordance with the requirements of 3.1.1.1 and 3.1.1.2.</p>				<p>E-1.4 Quality Control Verification</p> <p>The passivation procedure shall be performed in accordance with the requirements of 3.1.1.1 and 3.1.1.2. The passivation procedure shall be performed in accordance with the requirements of 3.1.1.1 and 3.1.1.2. The passivation procedure shall be performed in accordance with the requirements of 3.1.1.1 and 3.1.1.2.</p>				<p>E-1.5 Evaluation of Cleaned and Passivated Surfaces</p> <p>The passivation procedure shall be performed in accordance with the requirements of 3.1.1.1 and 3.1.1.2. The passivation procedure shall be performed in accordance with the requirements of 3.1.1.1 and 3.1.1.2. The passivation procedure shall be performed in accordance with the requirements of 3.1.1.1 and 3.1.1.2.</p>			

7. Rouging

- ISPE, ASME BPE, FDA, & USP

ISPE - whether or not the presence of rouge may be detrimental to the drug products. Refers to FDA, USP and ASME BPE

FDA - surfaces shall not be reactive, additive, or absorptive so as to alter the safety, identity, strength, quality, or purity of the drug product beyond the official or other established requirements.

USP - does not address design or material criteria directly, but rather indirectly by defining limits for the components that ultimately will enter the human body.

USP 788 & 789 Particulate matter in injections and ophthalmic solutions

	≥ 10 µm	≥ 25 µm
788 Injections	25 per mL	3 per mL
789 Ophthalmic	50 per mL	5 per mL

7 Rouging

– ASME BPE Nonmandatory appendix D, Rouge and stainless steel

- D-1 General
- D-2 Considerations for reducing rouge formation
 - Alloy selection/composition
 - Polishing
 - Passivation
 - Welding
 - Product form and fabrication
- D-3 Evaluation methods to measure rouge
- D-4 Methods to remediate the presence of rouge in a system

Table D-2-1 Considerations That Affect the Amount of Rouge Formation During the Fabrication of a System

Variables	Considerations
Category 3 — Strong Influence on the Formation of Rouge [Note (1)]	
Alloy selection	Selection of the proper alloy (e.g., 316L-type or 6 moly-type stainless steel) should address the corrosive effects of the process conditions. For example, low-carbon stainless steel (316L-type) has better corrosion resistance than higher-carbon stainless steels (316-type). Beneficial alloys can mitigate premature or accelerated corrosion. Higher nickel content will enhance corrosion resistance.
Mechanical polishing/buffing	Striations from cold working techniques may include residual grinding/polishing debris in lapping inclusions. Cumulative increase of interior area due to surface finish inconsistency proportionally exposes more alloy to the mechanisms of corrosion and burden of passivation.
Electropolishing	Minimizes the exposure area of the alloy to oxidizing fluids or halides and minimizes the origins for micropitting by corrosion mechanisms.
Passivation	Impedes or retards corrosive development of stainless steel surfaces. The effectiveness of passivation methods in terms of depth and enhancement of surface alloy ratios (i.e., chrome to iron) determines the eventual propensity of an alloy to corrode and the rate of corrosion.
Alloy composition	(% molybdenum, chromium, nickel, etc.) The microstructure quality affects precipitation of impurities at grain boundaries. Migration of impurities to the alloy surface can either support corrosion cells or seed downstream corrosion. Weld joints on tubing and/or components with dissimilar sulfur concentrations may result in lack of penetration due to weld pool shift. The resulting crevice may become a corrosion initiation site.
Welding, welding conditions, purging, etc.	Improper welds can result in chromium-depleted heat-affected zones (HAZs) and other conditions that reduce corrosion resistance. Weld discontinuities create opportunities to trap fluid-borne impurities. Cracks resulting from poor welds will create breaches in the passive layer and form active corrosion cells. Proper purging prevents weld contamination by heat tint oxides and the concurrent loss of corrosion resistance. Passivation cannot reverse the effects of improper purging.
Product form and fabrication methods	The ferrite content can be greatly affected by the forming process (e.g., a forging will typically have much lower ferrite percentages than a casting). Barstock endgrain voids at the surface can enhance the potential of the alloy to pit and corrode. Minimization of differences in sulfur content will enhance the potential for successful welding.
Category 2 — Moderate Influence on the Formation of Rouge [Note (1)]	
Installation/storage environment	Unidentified corrosion due to the storage or installation environment, including carbon steel contamination, scratching, exposure to chemical contaminants, stagnated condensation or liquids, etc., may warrant a derouging step prior to passivation. Failure to detect instances of corrosion will marginalize the effect of a normal passivation.
Expansion and modifications to an established system	Oxide formations in newly commissioned systems can form at different rates than in older systems and initially generate migratory Class I rouge. Where oxide films exist in established systems, they are likely to be more stable, producing less migratory iron or chrome oxides. Because the newer system can generate and distribute lightly held Class I migratory hematite forms throughout the system, the corrosion origin and cause can be difficult to identify.

NOTE: (1) There is well-established industry data supporting this, and it needs to be considered.

8. Installation & operation

- ISPE, FDA & ASME BPE

- Slope
- Space for inspection and service
- Pressure relief devices
- Over pressure
- Thermal shocks
- No drain of utility side


Table SD-2.4.3.1-1 Slope Designations for Gravity-Drained Lines

Slope Designation	Minimum Slope, in./ft	Minimum Slope, mm/m	Minimum Slope, %	Minimum Slope, deg
GSD1	$\frac{1}{16}$	5	0.5	0.29
GSD2	$\frac{1}{8}$	10	1.0	0.57
GSD3	$\frac{1}{4}$	20	2.0	1.15
GSD0	Line slope not required			

Ambient or cold system



Hot system



9. Service & maintenance

- ISPE & FDA



- Periodic disassembly and inspection
- Pressure testing
- Documentation
- Original spare parts

§ 211.67 Equipment cleaning and maintenance.

- (a) Equipment and utensils shall be cleaned, maintained, and, as appropriate for the nature of the drug, sanitized and/or sterilized at appropriate intervals to prevent malfunctions or contamination that would alter the safety, identity, strength, quality, or purity of the drug product beyond the official or other established requirements.
- (b) Written procedures shall be established and followed for cleaning and maintenance of equipment, including utensils, used in the manufacture, processing, packing, or holding of a drug product. These procedures shall include, but are not necessarily limited to, the following:
 - (1) Assignment of responsibility for cleaning and maintaining equipment;
 - (2) Maintenance and cleaning schedules, including, where appropriate, sanitizing schedules;
 - (3) A description in sufficient detail of the methods, equipment, and materials used in cleaning and maintenance operations, and the methods of disassembling and reassembling equipment as necessary to assure proper cleaning and maintenance;
 - (4) Removal or obliteration of previous batch identification;
 - (5) Protection of clean equipment from contamination prior to use;
 - (6) Inspection of equipment for cleanliness immediately before use.
- (c) Records shall be kept of maintenance, cleaning, sanitizing, and inspection as specified in §§ 211.180 and 211.182.

[43 FR 45077, Sept. 29, 1978, as amended at 73 FR 51931, Sept. 8, 2008]

Part 2: Guidance fulfilment and challenge solving

- With Alfa Laval Pharma-line S&T and Pharma-line Point of Use





Part 2: Guidance fulfilment and challenge solving

- With Alfa Laval Pharma-line S&T



Cross contamination and sanitization

- Double tubesheet design
- Tube expansion
- U-bend design
- Weld
- Weep hole / Drain slot
- Slope



U-bend



(f) Minimum recommended bend radii for heat exchangers should be as follows:

Nominal Tube O.D.		Minimum Bend Radius	
in.	mm	in.	mm
0.375	9.5	0.625	15.2
0.500	12.7	0.750	19.1
0.625	15.8	0.938	23.8
0.750	19.1	1.125	28.6
1.000	25.4	1.500	38.1

Nominal tube O.D.	Minimum Bend Radius
8	22
10	23
14	36

In Alfa Laval Pharma-line we are going beyond ASME BPE recommendation for bending radius



U-bend




Bad bending makes **SCC**
(Stress Corrosion Cracking)



In Alfa Laval Pharma-line we are going beyond the ASME BPE recommendation for bending radius

MJ-8.1 The weld shall not have any discontinuities such as cracks, voids, porosity, or joint misalignment that will promote contamination of the process fluid.

- Careful welding
- No grinding



- Grinding to correct bad welds causes weak construction and porosity
- Bad welding causes roughing





Weep hole / Drain slot

ASME BPE 2019 Chapter 2 : Design

- SD-3.6.1 G : Construction


Pass rib drain slot in bonnet




Slope




Pharma-line S →




Pharma-line P →




Slope



Pharma-line S →



Pharma-line P →



Inlet slope

Pharma-line S&T are sized for turbulent flow

HTRI	Tubeside Incremental Monitor							Page 1
	Released to the following HTRI Member Company: <i>Alfa Laval</i> <i>Alfa Laval</i>							
Xist 8.0.1 2021-11-04 15:30 SN: 00571-500098257							SI Units	
9615012521 BEU14-16-3								
Rating - Horizontal Multipass Flow Small BEU Shell With Single-Segmental Baffles								
Point number	(--)	1	2	3	4	5	6	
Tube Pass	(--)	1	1	1	1	1	1	
Local Reynolds	(--)	26815	26826	25608	24722	23881	23082	
Vapor Reynolds	(--)	0,0000	0,0000	0,0000	0,0000	0,0000	0,0000	
Liquid Reynolds	(--)	27614	26826	25608	24722	23881	23082	

		Page 6
		SI Units
	35	36
	2	2
	11673	11668
	0,0000	0,0000
	11673	11543

Surface finish



Surface Measurement Report of all different pipe for Pharma-line

Surface Ra-measurement pipe before bending and after bending.

- We have cut 2 pipe 8x1, 3 pipe 10x1 and 3 pipe 14x1 and taken RA- measurement for all pipes. From the beginning all pipes are RA<0,4µm Electropolished
- We have taken RA- measurement before and after the bending. In the protocol below can you see radius of the bending.
- After the bending all the pipes are within the specified surface finish.



Example

Product	Before bending 1 point	After bending 3 point
Rör 8x1 Radie 22 (inner)	0.119µm	0.376µm / 0.373µm / 0.145µm
Rör 8x1 Radie 32 (outer)	0.099µm	0.199µm / 0.085µm / 0.244µm
Rör 10x1 Radie 22 (inner)	0.132µm	0.210µm / 0.165µm / 0.330µm
Rör 10x1 Radie 34 (middle)	0.104µm	0.266µm / 0.176µm / 0.311µm
Rör 10x1 Radie 46 (outer)	0.074µm	0.152µm / 0.331µm / 0.242µm
Rör 14x1 Radie 36 (inner)	0.163µm	0.222µm / 0.149µm / 0.231µm
Rör 14x1 Radie 52 (middle)	0.234µm	0.260µm / 0.235µm / 0.267µm
Rör 14x1 Radie 67 (ytter)	0.081µm	0.292µm / 0.160µm / 0.183µm

Date: 2013-10-25

Classification of signature

Kristina Olsson
Quality Manager
Zetterströms Rostfria AB

Table SF-2.2-1 Acceptance Criteria for Metallic Process Contact Surface Finishes

Anomaly or Indication	Acceptance Criteria
Pits/porosity	If diameter <0.020 in. (0.51 mm) and bottom is shiny [Notes (1) and (2)], Pits <0.003 in. (0.08 mm) diameter are irrelevant and acceptable.
Cluster of pits/porosity	No more than 4 pits per 0.5 in. (13 mm) x 0.5 in. (13 mm) inspection window. The cumulative total diameter of all relevant pits shall not exceed 0.040 in. (1.02 mm).
Dents	None accepted [Note (3)]
Finishing marks	If R_a max. is met
Welds	Welds used in the as-welded condition shall meet the requirements of MJ-8. Welds finished after welding shall be flush with the base metal, and concavity and convexity shall meet the requirements of MJ-8. Such finishing shall meet the R_a requirements of Table SF-2.4.1-1.
Nicks	None accepted
Scratches	For tubing, if cumulative length is <12.0 in. (305 mm) per 20 ft (6.1 m) tube length or prorated and if depth is <0.003 in. (0.08 mm) For fittings, valves, and other process components, if cumulative length is <0.25 in. (6.4 mm), depth <0.003 in. (0.08 mm), and R_a max. is met For vessels, if length <0.50 in. (13 mm) at 0.003 in. (0.08 mm) depth and if <3 per inspection window [Note (4)]
Surface cracks	None accepted
Surface inclusions	If R_a max. is met
Surface residuals	None accepted, visual inspection
Surface roughness (R_a)	See Table SF-2.4.1-1
Weld slag	For tubing, up to 3 per 20 ft (6.1 m) length or prorated, if <75% of the width of the weld bead For fittings, valves, vessels, and other process components, none accepted (as welded shall meet the requirements of MJ-8 and Table MJ-8.4-1)
Blistering	None accepted

GENERAL NOTE: This table covers surface finishes that are mechanically polished or any other finishing method that meets the R_a max.


- NOTES:
- (1) Black bottom pit of any depth is not acceptable.
 - (2) Pits in superaluminic and nickel alloys may exceed this value. Acceptance criteria for pit size shall be established by agreement between owner/user and supplier. All other pit criteria remain the same.
 - (3) For vessels, dents in the area covered by and resulting from welding dimple heat transfer jackets are acceptable.
 - (4) An inspection window is defined as an area 4 in. x 4 in. (100 mm x 100 mm).



Comprehensive documentation




-  Certificate of design-examination, approved drawings and calculations
-  List of welders + Weld procedure + Welder qualification
-  Material certificates for pressure and product referred parts, 3.1 Certificate
-  Component certificates
-  Liquid Penetrant Test Report and procedure
-  Surface treatment test report + Ra values (if demand)
-  Dimension control report
-  Pressure test certificate (signed by Notified body)
-  Identification (technical sign)
-  CE-Documents
-  Installation, Operation and maintenance manual
-  Quality System certificates (ISO9000:2000, ASME)



Alfa Laval support



- Selection & Sizing
- Installation
- Operation
- Service & Maintenance
- Inspection
- Cleaning
- Etc.



Part 2: Guidance fulfilment and challenge solving

- Pharma-line Point of Use



Cross contamination and sanitization

- No internal welds
- U-bend design
- Mounting vertical
- Drain in lower points
- No dead leg, $L/D < 2$
- High Reynold number > 4000
- Self-sanitization in standby mode





Figure 6 PoU200 and 2" pitot tube connected to test rig, TR1 for Stand-by-mode test.



Video

A close-up photograph of a hand placing a white puzzle piece into a larger puzzle piece. The larger piece is red and has the word 'ANSWER' written on it in white. The smaller piece is white and has the word 'QUESTION' written on it in black. The background consists of many other white puzzle pieces, some of which are slightly out of focus. The hand is positioned on the right side of the frame, with fingers gripping the edges of the 'QUESTION' piece as it is being inserted into the 'ANSWER' piece.

QUESTION

ANSWER

Contacts

- For questions or support



Uday Prabhakar, Senior Inside Sales Engineer, Alfa Laval Dubai

Mobile : +971 56 996 2256

uday.prabhakar@alfalaval.com

Adnan Raza, Sales and Partner Sales Manager, Alfa Laval Dubai

Mobile: +971 56 407 2797

adnan.raza@alfalaval.com

Agus R. Adisuwondo, Regional Pharma Manager - SEA, Alfa Laval Indonesia

Mobile: +62 8111 351 0301

agus.adisuwondo@alfalaval.com

מחנה
התעסוקה